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
www.sageair.com
© 2011 Sage Air Technologies, Inc. All rights reserved.

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Reprocessing of flexible endoscopes and endoscopic accessories used in gastrointestinal endoscopy: Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology Nurses and Associates (ESGENA) – Update 2018

GUIDELINE

ASGE guideline for infection control during GI endoscopy

Prepared by: ASGE QUALITY ASSURANCE IN ENDOSCOPY COMMITTEE

Audrey H. Calderwood, MD, Lukejohn W. Day, MD, V. Raman Muthusamy, MD, James Collins, RN, Ralph David Hambrick, III, RN, Andrew S. Brock, MD, Nalini M. Guda, MD, Jonathan M. Buscaglia, MD, Bret T. Petersen, MD, Navtej S. Buttar, MD, Lauren G. Khanna, MD, Vladimir M. Kushnir, MD, Aparna Repaka, MD, Nicolas A. Villa, MD, Glenn M. Eisen, MD, MPH, Chair

This document was reviewed and approved by the Governing Board of the American Society for Gastrointestinal Endoscopy (ASGE).

5.1 Principles of infection control

RECOMMENDATION

As the carrier status of patients is often unknown, all patients should be treated as potentially infectious.

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► **Table 1** Spaulding classification and reprocessing of medical devices.

Spaulding classification	Examples in GI endoscopy	Reprocessing
Noncritical devices	<ul style="list-style-type: none"> ▪ Fingertip for pulse oximetry ▪ Blood pressure cuff ▪ Electrodes for high frequency surgery and ECG ▪ Mouthguard 	<ul style="list-style-type: none"> ▪ Manual cleaning and disinfection attaining at least a given level of bactericidal and yeasticidal activity
Semicritical devices	<ul style="list-style-type: none"> ▪ Flexible endoscopes and their endoscope components 	<ul style="list-style-type: none"> ▪ Thorough manual cleaning including brushing is mandatory, followed by: Reprocessing, including cleaning, disinfection (attaining at least a given level of minimum bactericidal, fungicidal, mycobactericidal, and virucidal activity), and rinsing ▪ Automated reprocessing in an EWD is strongly recommended ▪ Thorough drying before storage in closed cabinets or storage cabinets with a drying function <p>Competent staff specially trained in endoscope reprocessing (in line with national laws and regulations) are required.</p>
Critical devices	<ul style="list-style-type: none"> ▪ Endoscopic accessories, e. g. biopsy forceps, polypectomy snares, ERCP accessories, etc. ▪ Flexible endoscopes only if medical indication for sterilization is given 	<p>For reusable devices, validated and standardized reprocessing, preferably in a CSSD is strongly recommended, including:</p> <ul style="list-style-type: none"> ▪ Thorough cleaning ▪ Automated reprocessing systems ▪ Sterile packages ▪ Sterilization <p>Proof of structured training for reprocessing medical devices (in line with national laws and regulations)</p>

GI, gastrointestinal; ECG, electrocardiogram; EWD, endoscope washer-disinfector; ERCP, endoscopic retrograde cholangiopancreatography; CSSD, central sterilization and supply department.

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



Precauciones universales

RECOMMENDATION

All endoscopes and reusable
should be reprocessed with a u
processing procedure following
cedure (universal precautions).

RECOMMENDATION

A traceability system should be
patients in the case of an outbre

RECOMMENDATION

The endoscopy department sh
the carrier status of the patient, so any pertinent precau-
tions can be taken.

RECOMMENDATION

reprocessing procedure should
l protective equipment (PPE)

single-use gloves (EN 374);
sses or visors), face masks, and
hair covering;

equipment (RPE) when handling
disinfectants containing respira-

-resistant protection gowns (EN

ed throughout the entire repro-
er to avoid contact with infec-
and disinfectants.

✓ SANGRE

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✓ PIEL NO INTACTA

✓ MUCOSAS


¡IMPORTANTE!

1

2

Precautions in the endoscopy unit

A number of essential precautions should be observed in the endoscopy unit in order to minimize infectious risks to both personnel and patients. Hands should be washed before and after each patient interaction, whether or not gloves are worn. The use of soap and water is required when hands are visibly soiled or an employee has an encounter with a patient with a suspected and/or known infectious cause of diarrhea. In all other cases, alcohol-based agents are acceptable.^{102,126} In endoscopy units, the prevention of *C difficile* transmission should be considered when endoscopy is performed on patients with diarrhea or known *C difficile* infection. Handwashing with soap and water should be undertaken for mechanical removal of spores from employee hands. Similarly, the

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RECOMMENDATION

Outbreaks should be managed within the multidisciplinary team of endoscopy departments, hospital hygiene experts, microbiologists, manufacturers, and regulatory bodies, if applicable.

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www.elsevier.es/eimc

Review article

Microbiological monitoring of medical devices after cleaning, disinfection and sterilisation[☆]

Rosa María Blázquez-Garrido^{a,*}, Eva Cuchí-Burgos^b, Carmen Martín-Salas^c, Patricia Ruiz-Garbajosa^d

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Table 1

Recommendations for microbiological study of endoscopes.

Cut-off point	Microorganisms	Frequency	Guidelines
<10 CFUs per device	Low-risk microbiota (skin and environment) Coagulase-negative <i>Staphylococcus</i> Diphtheroids <i>Micrococcus</i> spp. <i>Bacillus</i> spp.	Not established (except in duodenoscopes: every 60 procedures or monthly)	CDC ⁹
Any count	High-risk microbiota <i>S. aureus</i> <i>Enterococcus</i> spp. <i>S. viridans</i> <i>P. aeruginosa</i> Enterobacteria		
Same cut-off points as ESGENA and GESA	<i>Enterococcus</i> spp. Enterobacteria <i>S. viridans</i> NF GNB Mycobacteria on bronchoscopes Atypical NF GNB mycobacteria and <i>Legionella</i> spp. in water and on washers	Monthly; once circuit is established, may be quarterly	SEIMC 2012 ²⁰ UNE EN-ISO 15883 ²¹
No growth 1–10 CFUs 10–100 CFUs 100–1000 CFUs > 10,000 CFUs	<i>Enterococcus</i> spp. Enterobacteria <i>S. viridans</i> NF GNB Rapidly growing mycobacteria on bronchoscopes	Monthly Duodenoscopes Bronchial endoscopes Other flexible endoscopes Quarterly Other gastrointestinal endoscopes	GESA 2010 ¹⁸
<20 CFUs/channel Indicator microorganisms; any count	Any type of microbiota Indicator microorganisms Enterobacteria <i>P. aeruginosa</i> <i>S. aureus</i> Final wash water: Rapidly growing mycobacteria <i>Legionella</i> spp.	Quarterly or more often	ESGE-ESGENA 2007 ¹⁶
Microorganisms not included Anaerobes Viruses <i>Helicobacter pylori</i>		No guidelines mention fungi	

CDC: Centres for Disease Control and Prevention; CFUs: colony-forming units; ESGE-ESGENA: European Society of Gastrointestinal Endoscopy-European Society of Gastroenterology and Endoscopy Nurses and Associates; GESA: Gastroenterological Society of Australia; NF GNB: non-fermenting Gram-negative bacilli; SEIMC: Spanish Association of Infectious Diseases and Clinical Microbiology.

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Table 2

General interpretation of microbiological results.


Microorganisms isolated	Interpretation/action to be taken
Growth of <i>S. epidermidis</i> in low counts	Suspected contamination in sampling or sample processing or storage It is recommended that sampling be repeated
Significant counts of microorganisms on several devices Isolation of enteropathogens such as <i>Salmonella</i> spp. Isolation of <i>P. aeruginosa</i>	Suspected failures in the cleaning/disinfection of devices process It is recommended that the cleaning and disinfection procedure be reviewed and that sampling be repeated
Isolation of <i>P. aeruginosa</i> and other non-fermenting GNB	Suspected failures in drying or storage of devices or contamination of wash water. Assess the suitability of placing filters
Repeated growth in significant counts of a microorganism on a single device	Suspected structural problems with the devices It is recommended that the manufacturer review the device and that any leaks be managed

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RECOMMENDATION

The room should have:

- Appropriate size and lighting, and ventilation and fume extraction in order to minimize the risks from chemical vapors;
- Appropriate technical equipment and protective measures in order to ensure safe reprocessing following standardized and validated reprocessing procedures;
- Strict spatial or at least operational separation of dirty and clean/storage areas, in order to avoid recontamination of reprocessed endoscopes and endoscopic accessories.

This should be supported by the room architecture and design as well as by the one-way workflow from dirty to clean areas. Ideally, the standards should comply with those of the central sterilization and supply department (CSSD) in the particular country.

RECOMMENDATION

It is the responsibility of the clinical service provider to ensure that adequate facilities for reprocessing are available.

RECOMMENDATION

Independently of the distance between endoscopy rooms and reprocessing area, the workflow should ensure immediate reprocessing of used equipment.

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7.1 Documentation

RECOMMENDATION

The complete reprocessing cycle should be documented:

- Each reprocessing step (including bedside cleaning, manual cleaning, and automated reprocessing in an EWD or ADD) should be recorded manually or electronically, including the names of the persons undertaking each step.
- The process parameters of the EWD and storage cabinets should be documented by printouts or electronically.
- All endoscopes should have a record of their reprocessing showing that they are ready for use on patients.
- The reprocessing record should be documented in the patient's files.

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