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Reprocessing of flexible endoscopes and endoscopic accessories used in gastrointestinal endoscopy: Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology Nurses and Associates (ESGENA) – Update 2018

Authors

Ulrike Bellenhoff¹, Holger Biering², Reinhard Blum³, Jadranka Brjalic⁴, Monika Cimbro⁵, Jean-Marc Dumonceau⁶, Cesare Hassan⁷, Michael Jung⁸, Birgit Kamp⁹, Christiane Neumann⁹, Michael Pietsch¹⁰, Lionel Pinau¹¹, Thierry Ponchon¹², Stanislav Rajceta¹³, Jean-François Rey¹⁴, Verona Schmidt¹⁵, Jayne Tillet¹⁶, Jeanin E. van Hooft¹⁷

Institutions

- 1 ESCENA Scientific Secretary, Ulm, Germany
- 2 Grevenbroich, Germany
- 3 Olympus Europe, Hamburg, Germany
- 4 University Hospital IBC-Zagreb-Rebro, Zagreb, Croatia
- 5 Medical Devices Division, CIRC Europe Srl, Nova Milanese (MB), Italy
- 6 Gedyt Endoscopy Center, Buenos Aires, Argentina
- 7 Digestive Endoscopy Unit, Catholic University, Rome, Italy
- 8 2nd Department of Internal Medicine, Katholisches Klinikum, Mainz, Germany
- 9 ESCENA Past President, Birmingham, UK
- 10 Department of Hygiene and Infection Prevention, Medical Center, University Hospital, Mainz, Germany
- 11 Biotech Germande, Marseille, France
- 12 Digestive Diseases Department, Hôpital Edouard Herriot, Lyon, France
- 13 2nd Department of Internal Medicine, Charles University Teaching Hospital, Hradec Králové, Czech Republic
- 14 Institut Arnauld Tzanck, St. Laurent du Var, France
- 15 Microbiology and Hygiene Department, Chemische Fabrik Dr. Weigert, Hamburg, Germany
- 16 St. Woolos Hospital, Newport, UK
- 17 Department of Gastroenterology and Hepatology, Amsterdam University Medical Centers, University of Amsterdam, The Netherlands

Bibliography

DOI: <https://doi.org/10.1055/a-0799-1629>
 Published online: 20.11.2018 | Endoscopy 2018; 50: 1205–1234
 © Georg Thieme Verlag KG Stuttgart · New York
 ISSN 0013-726X

Corresponding author

Ulrike Bellenhoff, ESCENA Scientific Secretary, Ferdinand-Sauerbruch-Weg 16, 89075 Ulm, Germany
info@esgena.org

ABSTRACT

This Position Statement from the European Society of Gastrointestinal Endoscopy (ESGE) and the European Society of Gastroenterology Nurses and Associates (ESGENA) sets standards for the reprocessing of flexible endoscopes and endoscopic devices used in gastroenterology. An expert working group of gastroenterologists, endoscopy nurses, chemists, microbiologists, and industry representatives reprocessing in order to maintain hygiene and infection control.

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

 CANTEL MEDICAL Cantel Medical (Italy) S.r.l. Single sharador Via Laurentina, 169 00071 Pomezia (RM) - Italy Tel.: +39 06 9145399 Fax: +39 06 9146099	TEST CENTRE GLP Certified by the <i>Ministry of Health</i> according to DL 02/03/07 N° 50- Directive 2004/9/EC for Good Laboratory Practice	
	Study No. 05.15	Final Report of 07/01/2015 COPY IN ENGLISH N°  ON 10/12/2015

 CANTEL MEDICAL Cantel Medical (Italy) S.r.l. Single sharador Via Laurentina, 169 00071 Pomezia (RM) - Italy Tel.: +39 06 9145399 Fax: +39 06 9146099	TEST CENTRE GLP Certified by the <i>Ministry of Health</i> according to DL 02/03/07 N° 50- Directive 2004/9/EC for Good Laboratory Practice	
	Study No. 05.15	Final Report of 07/01/2015 COPY IN ENGLISH N°  ON 10/12/2015

VALIDATION OF THE STERILIZING PROCESS

IN THE WASHER- DISINFECTOR - STERILIZER

MOD. MEDIVATORS® ISA™ (S/ N° 6-0001)

ACCORDING TO STANDARD UNI EN ISO 14937:2009

CLIENT/PLACE OF TEST

CANTEL MEDICAL (ITALY) S.r.l. – via Laurentina, 169 – 00071 Pomezia

TEST PERIOD:

FROM: 07/29/15 to 07/31/2015

REFERENCE STANDARD

- **UNI EN ISO14937:2009:** "Sterilization of health care products- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices (ref. par 6.2, 6.3, 7)
- **Annex D: Approach 3- UNI EN ISO 14937:2009:** "Conservative process definition based on inactivation of reference microorganisms."
- **Ph EU: 8.0 01/2014- par 2.6.1:** "Biological tests: sterility"
- **UNI EN ISO 11138-1:2006:** Sterilization of health care products- Biological indicators- Part 1: General requirements.

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